









Ø0.13in [3.2mm]  
SLIP FIT HOLES FOR  
FIBERBENCH ACCESSORIES  
10 PLACES

1/4-20(M6) CLEARANCE SLOT  
FOR MOUNTING  
2 PLACES



ISOMETRIC VIEW  
SCALE: 3:4

FOR INFORMATION ONLY  
NOT FOR MANUFACTURING PURPOSES

|                                                                       |      |           |                                                                                                           |               |       |
|-----------------------------------------------------------------------|------|-----------|-----------------------------------------------------------------------------------------------------------|---------------|-------|
| DRAWING PROJECTION                                                    |      |           | <br>www.thorlabs.com |               |       |
|                                                                       | NAME | DATE      |                                                                                                           |               |       |
| DRAWN                                                                 | CB   | 08/AUG/11 | FIBER BENCH, 51mm                                                                                         |               |       |
| APPROVAL                                                              | BAG  | 12/SEP/18 | MATERIAL                                                                                                  |               | REV E |
| COPYRIGHT © 2018 BY THORLABS                                          |      |           | N/A                                                                                                       |               |       |
| VALUES IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUND OFF ERRORS |      |           | ITEM #                                                                                                    | APPROX WEIGHT |       |
|                                                                       |      |           | FB-51                                                                                                     | 0.22 kg       |       |