




ISOMETRIC VIEW
SCALE 1:1

SM1 (1.035"-40) SERIES
INTERNAL THREAD
SM1RR RETAINING RING INCLUDED
CLEAR APERTURE 0.09in [22.9mm]



INCLUDED BUT NOT SHOWN:

- FOUR #2-56 SOCKET HEAD SCREW
- TWO #8-32 SOCKET HEAD SHOULDER SCREW
- ONE SHORT ARM HEX-L-KEY 5/64
- ONE SHORT ARM HEX-L-KEY 3/32

FOR INFORMATION ONLY
NOT FOR MANUFACTURING PURPOSES

DRAWING PROJECTION					THORLABS www.thorlabs.com	
	NAME	DATE	FIBERBENCH WALL PLATE, SM1 THREAD			
DRAWN	JJM	18/JUL/08	MATERIAL			
APPROVAL	KAH	13/FEB/13	303 STAINLESS STEEL			REV C
COPYRIGHT © 2012 BY THORLABS			ITEM #			
VALUES IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUND OFF ERRORS			HCA3-SM1		APPROX WEIGHT 0.06 kg	