




Item #	Color	Correlated Color Temperature	Minimum Power LED Output *	Typical Power LED Output *	Test Current for Typical LED Power	Maximum Current (CW)	Forward Voltage	Bandwidth (FWHM)	Typical Lifetime
MCWHF2	Cold White	6200 K	21.5 mW	27.0 mW	1000 mA	1000 mA	2.9 V	N/A	>50 000 h

* for MM Fiber with 400 μ m core diameter, 0.39 NA (FT400EMT)

FOR INFORMATION ONLY
NOT FOR MANUFACTURING PURPOSES

DRAWING PROJECTION



	NAME	DATE
DRAWN	Hit.	24/FEB/16
APPROVAL	KB	24/FEB/16

COPYRIGHT © 2016 BY THORLABS

VALUES IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUND OFF ERRORS

THORLABS
www.thorlabs.com

COLD WHITE (6200 K)
FIBER-COUPLED LED

MATERIAL

N/A

REV
A

ITEM #

MCWHF2

APPROX WEIGHT

0.12 kg