






FOR INFORMATION ONLY
NOT FOR MANUFACTURING PURPOSES

DRAWING PROJECTION					THORLABS www.thorlabs.com	
	NAME	DATE				
DRAWN	WZ	29/MAR/18	$\varnothing 6.0\text{mm}$ MOUNTING POST 25.0mm LONG			REV C
APPROVAL	AK	11/APR/18	MATERIAL			
COPYRIGHT © 2018 BY THORLABS			STAINLESS STEEL			
VALUES IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUND OFF ERRORS			ITEM #		APPROX WEIGHT	
			MS1R/M		5.3 g	