



REMOVABLE SS4S025 (INCLUDED)
(4-40 SETSCREW, 0.25in [6.35mm] LONG)
0.050in [1.3mm] HEX
2 PLACES



FOR INFORMATION ONLY
NOT FOR MANUFACTURING PURPOSES

DRAWING PROJECTION				THORLABS www.thorlabs.com	
	NAME	DATE	MINI SERIES CAGE ROD, 4.0" LONG PACK OF 4		
DRAWN	KBH	11/OCT/16			
APPROVAL	AK	12/OCT/16			
COPYRIGHT © 2016 BY THORLABS			MATERIAL		REV
			STAINLESS STEEL		A
VALUES IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUND OFF ERRORS			ITEM #		APPROX WEIGHT
			SR4-P4		0.01 kg