DDSA01 (/M)

End Mounting Plate for DDS220 Stages



Original Instructions

Contents

Chaper 1 Overview 1
1.1 Description *1*
Chaper 2 Safety 2
2.1 Safety Information *2*
Chaper 3 Installation 3
3.1 Installation *3*
3.2 Typical Application *3*
Chaper 4 Specifications 4
Chaper 5 Thorlabs Worldwide Contacts 5

Chapter 1 Overview

1.1 Description

The DDSA01 adapter plate fixes to either end of the DDS220 stage and allows accessories and device holders to be mounted. The plates are drilled with an array of 1/4”-20 (M6) holes and the finish is a black, low reflective anodized coating.

Chapter 2 Safety

2.1 Safety Information

For the continuing safety of the operators of this equipment, and the protection of the equipment itself, the operator should take note of the **Warnings, Cautions** and **Notes** throughout this handbook and, where visible, on the product itself.

The following safety symbols may be used throughout the handbook and on the equipment itself.

The following safety symbols may be used throughout the handbook and on the equipment itself.



Warning: Risk of Electrical Shock

Given when there is a risk of injury from electrical shock.



Warning

Given when there is a risk of injury to users.



Caution

Given when there is a risk of damage to the product.

Note

Clarification of an instruction or additional information.

Chapter 3 Installation

3.1 Installation




Fig. 3.1 Installation

Fig. 3.2 Installation:

1) Using two 1/4”-20 x 3/8 (M6 x 16) bolts (not supplied), fit the DDSA01 plate to either end of the DDS220 stage.

3.2 Typical Application

When mounted at either end of the DDS220 stage, the mounting plates make an ideal platform for mounting accessories for applications such as a delay line.



Fig. 3.3 DDS220 stage with DDSA01 plates fitted

Chapter 4 Specifications

Specification	Value
Working Height of Stage Assembly	44.4 mm (1.75”)
Mounting Area Dimensions	90 mm x 64.0 mm x 24.0 mm (3.54" x 2.52" x 0.94”)
Weight:	0.33 kg (0.72 lb)
Construction	Solid Aluminum
Finish	Black Anodized

Chapter 5 Thorlabs Worldwide Contacts

For technical support or sales inquiries, please visit us at www.thorlabs.com/contact for our most up-to-date contact information.



USA, Canada, and South America
Thorlabs, Inc.
sales@thorlabs.com
techsupport@thorlabs.com

Europe
Thorlabs GmbH
europe@thorlabs.com

France
Thorlabs SAS
sales.fr@thorlabs.com


Japan
Thorlabs Japan, Inc.
sales@thorlabs.jp

UK and Ireland
Thorlabs Ltd.
sales.uk@thorlabs.com
techsupport.uk@thorlabs.com

Scandinavia
Thorlabs Sweden AB
scandinavia@thorlabs.com

Brazil
Thorlabs Vendas de Fotônicos Ltda.
brasil@thorlabs.com

China
Thorlabs China
chinasales@thorlabs.com



THORLABS

www.thorlabs.com